

## IPG Cantech Process Indicator Tape 144C

### EXPOSURE INTERPRETATION:

Steam Sterilization Type 1

Tape is a single use device that is interpreted for appropriate signal color after processing.

### Before Cycle: Unexposed 144C



*Designed for use in gravity and pre-vacuum sterilizers*

### Cycle Completed: Exposed 144C



*Indicator lines on the tape turn black when exposed to proper steam sterilization conditions.*

*Lead is not used in the manufacturing of the ink of indicator lines.*



*Other examples that show acceptable "exposed" results*